

À des fins de recherche uniquement

Anticorps Monoclonal anti-RAB3D

Numéro de catalogue: 67118-1-Ig



Informations de base

Numéro de catalogue: 67118-1-Ig	Numéro d'acquisition GenBank: BC016471	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 9545	CloneNo.: 3C2F10
Hôte: Mouse	Nom complet: RAB3D, member RAS oncogene family	Dilutions recommandées: WB 1:5000-1:20000 IHC 1:500-1:2000 IF 1:400-1:1600
Isotype: IgG2a	MW calculé 219 aa, 24 kDa	
Immunogen Catalog Number: AG16784	MW observés: 28 kDa	

Applications

Applications testées: IF, IHC, WB, ELISA	Contrôles positifs:
Spécificité de l'espèce: Humain	WB : cellules NCI-H1299, cellules A549, cellules HT-29, cellules THP-1
Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.	IHC : tissu de cancer du côlon humain, IF : cellules A549,

Informations générales

Stockage


Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




NCI-H1299 cells were subjected to SDS PAGE followed by western blot with 67118-1-Ig (RAB3D antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67118-1-Ig (RAB3D antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67118-1-Ig (RAB3D antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using RAB3D antibody (67118-1-Ig, Clone: 3C2F10) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).