

À des fins de recherche uniquement

Anticorps Monoclonal anti-RASEF

Numéro de catalogue: 67159-1-Ig



Informations de base

Numéro de catalogue: 67159-1-Ig	Numéro d'acquisition GenBank: BC023566	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 2000 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 158158	CloneNo.: 1C2B2
Hôte: Mouse	Nom complet: RAS and EF-hand domain containing	Dilutions recommandées: WB 1:1000-1:6000 IF 1:200-1:800
Isotype: IgG2a	MW calculé 740 aa, 83 kDa	
Immunogen Catalog Number: AG2149	MW observés: 90 kDa	

Applications

Applications testées: IF, WB, ELISA	Contrôles positifs:
Spécificité de l'espèce: Humain	WB : cellules LNCaP, IF : cellules A431,

Informations générales

The RASEF (RAS and EF-hand domain-containing protein), a member of Rab GTPase protein family, contains a Rab GTPase domain in the C-terminal region, a coiled coil motif in an internal region as well as 2 EF hand domains which were of importance for binding to calcium ions in the N-terminus.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



LNCaP cells were subjected to SDS PAGE followed by western blot with 67159-1-Ig (RASEF antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using RASEF antibody (67159-1-Ig, Clone: 1C2B2) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).