

À des fins de recherche uniquement

Anticorps Monoclonal anti-ABTB1

Numéro de catalogue: 67246-1-Ig



Informations de base

Numéro de catalogue: 67246-1-Ig	Numéro d'acquisition GenBank: BC011858	Méthode de purification: Purification par protéine G
Taille: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 80325	CloneNo.: 2D2B11
Hôte: Mouse	Nom complet: ankyrin repeat and BTB (POZ) domain WB 1:1000-1:4000 containing 1	Dilutions recommandées:
Isotype: IgG1	MW calculé 54 kDa	
Immunogen Catalog Number: AG28959		

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : tissu de cervelet de porc,
Spécificité de l'espèce: Humain, porc	

Informations générales

Stockage


Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



pig cerebellum tissue were subjected to SDS PAGE followed by western blot with 67246-1-Ig (ABTB1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.