

À des fins de recherche uniquement

Anticorps Monoclonal anti-HOMER2

Numéro de catalogue: 67377-1-Ig



Informations de base

Numéro de catalogue: 67377-1-Ig	Numéro d'acquisition GenBank: BC012109	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1600 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 9455	CloneNo.: 1G5G1
Hôte: Mouse	Nom complet: homer homolog 2 (Drosophila)	Dilutions recommandées: WB 1:1000-1:6000 IHC 1:200-1:800
Isotype: IgG3	MW calculé 41 kDa	
Immunogen Catalog Number: AG25177	MW observés: 40 kDa	

Applications

Applications testées: IHC, WB, ELISA	Contrôles positifs:
Spécificité de l'espèce: Humain	WB : cellules LNCaP, cellules HEK-293 IHC : tissu de cancer de l'endomètre humain,
Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6.0.	

Informations générales

Stockage


Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



LNCaP cells were subjected to SDS PAGE followed by western blot with 67377-1-Ig (HOMER2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissue slide using 67377-1-Ig (HOMER2 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissue slide using 67377-1-Ig (HOMER2 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).