

À des fins de recherche uniquement

# Anticorps Monoclonal anti-PIP5K1A

Numéro de catalogue: 67483-1-Ig



## Informations de base

Numéro de catalogue: 67483-1-Ig	Numéro d'acquisition GenBank: BC007833	Méthode de purification: Purification par protéine A
Taille: 150ul , Concentration: 2000 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 8394	CloneNo.: 3E5G12
Hôte: Mouse	Nom complet: phosphatidylinositol-4-phosphate 5- kinase, type I, alpha	Dilutions recommandées: WB 1:5000-1:50000 IHC 1:250-1:1000
Isotype: IgG2b	MW calculé 56 kDa, 63 kDa	
Immunogen Catalog Number: AG8355	MW observés: 63 kDa	

## Applications

### Applications testées:

IHC, WB, ELISA

### Spécificité de l'espèce:

Humain

**Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (\*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.**

### Contrôles positifs:

WB : cellules LNCaP, cellule HeLa, cellules HEK-293, cellules HeLa, cellules HepG2, cellules Jurkat, cellules K-562, cellules L02, cellules MOLT-4

IHC : tissu d'amygdalite humain,

## Informations générales

Phosphatidylinositol 4-phosphate 5-kinases (PIP5Ks) play diverse roles in the cellular biology of many organisms, including signal transduction, secretion and vesicular trafficking, and regulation of cytoskeleton assembly (PMID:17688436). There are three PIP5K isoforms,  $\alpha$ ,  $\beta$ , and  $\gamma$ . The nomenclature for the  $\alpha$  and  $\beta$  isoforms is switched between humans and mice (PMID:22096541). There are also several splicing variants of the  $\gamma$  isozyme have been identified (PMID:20945365). PIP5K1A has some isoforms produced by alternative splicing with the molecular mass of 56-63 kDa.

## Stockage

### Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

### Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquote n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Données de validation sélectionnées




Various lysates were subjected to SDS PAGE followed by western blot with 67483-1-Ig (PIP5K1A antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Beta Actin Monoclonal antibody (HRP-66009) as loading control.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67483-1-Ig (PIP5K1A antibody) at dilution of 1:500 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 67483-1-Ig (PIP5K1A antibody) at dilution of 1:500 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Various lysates were subjected to SDS PAGE followed by western blot with 67483-1-Ig (PIP5K1A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.