

À des fins de recherche uniquement

Anticorps Monoclonal anti-AgBR1

Numéro de catalogue: 67743-1-Ig



Informations de base

Numéro de catalogue: 67743-1-Ig	Numéro d'acquisition GenBank: XM_001844620	Méthode de purification: Purification par protéine G
Taille: 150ul , Concentration: 1000 µg/ml by Nanodrop and 439 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 6034353	CloneNo.: 2D1F10
Hôte: Mouse	Nom complet: AgBR1	Dilutions recommandées: WB 1:20000-1:100000
Isotype: IgG1	MW calculé 49 kDa	
Immunogen Catalog Number: AG30168	MW observés: 40 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : Ag30168,
Spécificité de l'espèce: Moustique	

Informations générales

AgBR1, a mosquito salivary protein also known as A. aegypti bacteria-responsive protein 1, is found to induce inflammatory responses at the bite site (PMID: 32218189). Various studies suggest that AgBR1 plays a role in transmitting or facilitating infections of several viruses via Aedes aegypti mosquitoes (PMID: 30858571, 31312526).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 67743-1-Ig at dilution of 1:160000 incubated at room temperature for 1.5 hours.