

À des fins de recherche uniquement

Anticorps Monoclonal anti-DAXX

Numéro de catalogue: 67879-1-Ig



Informations de base

Numéro de catalogue: 67879-1-Ig	Numéro d'acquisition GenBank: BC109074	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 1616	CloneNo.: 1F3G9
Hôte: Mouse	Nom complet: death-domain associated protein	Dilutions recommandées: WB 1:5000-1:50000 IHC 1:500-1:2000 IF 1:200-1:800
Isotype: IgG2b	MW calculé: 740 aa, 81 kDa	
Immunogen Catalog Number: AG14606	MW observés: 100-130 kDa	

Applications

Applications testées:

IF, IHC, WB, ELISA

Spécificité de l'espèce:

Humain

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (*) À défaut, 'le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.

Contrôles positifs:

WB : cellules Raji, cellules Daudi, cellules HeLa, cellules Jurkat, cellules MOLT-4

IHC : tissu de cancer du sein humain,

IF : cellules A431,

Informations générales

Death domain-associated protein (DAXX) is a transcription repressor involved in both physiological and pathological conditions. DAXX overexpression is a common feature in diverse cancers, which correlates with tumorigenesis, disease progression and treatment resistance. Daxx migrated with an apparent molecular weight of approximately 100-130 kDa. (PMID: 31350900, PMID: 32203224, PMID: 31614769)

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



Raji cells were subjected to SDS PAGE followed by western blot with 67879-1-Ig (DAXX antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 67879-1-Ig (DAXX antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed A431 cells using DAXX antibody (67879-1-Ig, Clone: 1F3G9) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).