

À des fins de recherche uniquement

# Anticorps Monoclonal anti-PI3 Kinase p110 Delta



Numéro de catalogue: 67964-1-Ig

## Informations de base

|   |  |   |
|---|--|---|
| Numéro de catalogue:<br>67964-1-Ig                          | Numéro d'acquisition GenBank:<br>BC132919                                  | Méthode de purification:<br>Purification par protéine G         |
| Taille:<br>150ul, Concentration: 1000 µg/ml by<br>Nanodrop; | Identification du gène (NCBI):<br>5293                                     | CloneNo.:<br>2G11C6   |
| Hôte:<br>Mouse  | Nom complet:<br>phosphoinositide-3-kinase, catalytic,<br>delta polypeptide | Dilutions recommandées:<br>WB 1:2000-1:10000<br>IF 1:400-1:1600 |
| Isotype:<br>IgG1  | MW calculé<br>1044 aa, 120 kDa   |   |
| Immunogen Catalog Number:<br>AG26995                        | MW observés:<br>110 kDa  |   |

## Applications

|  |   |
|--|---|
| Applications testées:<br>IF, WB, ELISA | Contrôles positifs:<br>WB : cellules LNCaP, cellules Jurkat<br>IF : cellules LNCaP, |
| Spécificité de l'espèce:<br>Humain     |   |

## Informations générales

### Stockage

**Stockage:**  
Stocker à -20°C. Stable pendant un an après l'expédition.  
**Tampon de stockage:**  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Données de validation sélectionnées



LNCaP cells were subjected to SDS PAGE followed by western blot with 67964-1-Ig (PIK3CD antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Jurkat cells were subjected to SDS PAGE followed by western blot with 67964-1-Ig (PIK3CD antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed LNCaP cells using PIK3CD antibody (67964-1-Ig, Clone: 2G11C6) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).