

À des fins de recherche uniquement

# Anticorps Monoclonal anti-Flotillin 1

Numéro de catalogue: 67968-1-Ig



## Informations de base

|                           |  |                                |             |                          |                             |
|---------------------------|--|--------------------------------|-------------|--------------------------|-----------------------------|
| Numéro de catalogue:      | 67968-1-Ig                                     | Numéro d'acquisition GenBank:  | BC001146    | Méthode de purification: | Purification par protéine A |
| Taille:                   | 150ul , Concentration: 1000 µg/ml by Nanodrop; | Identification du gène (NCBI): | 10211       | CloneNo.:                | 1D6C6                       |
| Hôte:                     | Mouse  | Nom complet:                   | flotillin 1 | Dilutions recommandées:  | WB 1:5000-1:50000           |
| Isotype:                  | IgG2b  | MW calculé:                    | 47 kDa      |                          |                             |
| Immunogen Catalog Number: | AG8694   | MW observés:                   | 47 kDa      |                          |                             |

## Applications

### Applications testées:

WB, ELISA

### Spécificité de l'espèce:

Humain, Lapin, porc, rat, souris

### Contrôles positifs:

WB : cellules LNCaP, cellules Jurkat, cellules MOLT-4, cellules Raji, plaquettes de sang périphérique humain, tissu cérébral de lapin, tissu cérébral de porc, tissu cérébral de rat, tissu cérébral de souris

## Informations générales

Flotillins are membrane-associated proteins that contain two homologous isoforms, flotillin 1 (FLOT1) and flotillin 2 (FLOT2). They are highly conserved and widely expressed proteins that commonly used as markers of lipid rafts. Flotillins are involved in various signaling pathways, cell adhesion, membrane trafficking and axonal growth.

## Stockage

### Stockage:

Stocker à -20 °C.

### Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 67968-1-Ig (Flotillin 1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.