

À des fins de recherche uniquement

Anticorps Monoclonal anti-ACTR3B

Numéro de catalogue: 68037-1-Ig



Informations de base

Numéro de catalogue: 68037-1-Ig	Numéro d'acquisition GenBank: BC008682	Méthode de purification: Purification par protéine G
Taille: 150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 57180	CloneNo.: 1G12C6
Hôte: Mouse	Nom complet: ARP3 actin-related protein 3 homolog B (yeast)	Dilutions recommandées:
Isotype: IgG1	MW calculé 418 aa, 48 kDa	
Immunogen Catalog Number: AG10218	MW observés: 48 kDa	

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, Lapin, porc, rat, souris

Contrôles positifs:

WB : cellules MOLT-4, cellules Jurkat, cellules K-562, cerveau de lapin, cerveau de rat, cerveau de souris, tissu cérébral de porc

Informations générales

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various cells were subjected to SDS PAGE followed by western blot with 68037-1-Ig (ACTR3B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.



Brain tissues were subjected to SDS PAGE followed by western blot with 68037-1-Ig (ACTR3B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.