

À des fins de recherche uniquement

Anticorps Monoclonal anti-NBR1

Numéro de catalogue: 68062-1-Ig



Informations de base

Numéro de catalogue: 68062-1-Ig	Numéro d'acquisition GenBank: BC009808	Méthode de purification: Purification par protéine G
Taille: 150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 4077	CloneNo.: 1F2B2
Hôte: Mouse	Nom complet: neighbor of BRCA1 gene 1	Dilutions recommandées: WB 1:2000-1:10000
Isotype: IgG1	MW calculé: 966 aa, 107 kDa	
Immunogen Catalog Number: AG8652	MW observés: 140 kDa	

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Contrôles positifs:

WB : cellules NCCIT, cellules HEK-293, cellules HeLa, cellules hTERT-RPE1, cellules K-562, tissu testiculaire de rat, tissu testiculaire de souris

Informations générales

NBR1, also named as 1A13B, KIAA0049 and M17S2, acts probably as a receptor for selective autophagosomal degradation of ubiquitinated targets. NBR1 and P62 can bind to autophagic effector proteins (Atg8 in yeast, MAP1LC3 protein family in mammals) anchored in the membrane of autophagosomes. It is a highly conserved multidomain scaffold protein with proposed roles in endocytic trafficking and selective autophagy. NBR1 is a novel PB1 adapter in Th2 differentiation and asthma. It functions as an autophagy receptor involved in targeting ubiquitinated proteins for degradation. It also has a dual role as a scaffold protein to regulate growth-factor receptor and downstream signaling pathways. Observed MW of NBR1 is 140 kDa (PMID: 22654911, PMID: 22484440).

Stockage

Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68062-1-Ig (NBR1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.