

À des fins de recherche uniquement

Anticorps Monoclonal anti-RCL1

Numéro de catalogue: 68221-1-Ig



Informations de base

Numéro de catalogue:	68221-1-Ig	Numéro d'acquisition GenBank:	BC001025	Méthode de purification:	Purification par protéine G
Taille:	150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	10171	CloneNo.:	1D1B3
Hôte:	Mouse	Nom complet:	RNA terminal phosphate cyclase-like 1	Dilutions recommandées:	WB 1:5000-1:50000
Isotype:	IgG1	MW calculé	41 kDa		
Immunogen Catalog Number:	AG7058	MW observés:	41 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules MCF-7, cellules HEK-293, cellules HeLa, cellules HepG2, cellules HSC-T6, cellules Jurkat, cellules K-562, cellules NIH/3T3
Spécificité de l'espèce:	Humain, rat, souris		

Informations générales

RNA 3'-terminal phosphate cyclase-like protein (RCL1) does not have cyclase activity. It plays a role in 40S-ribosomal-subunit biogenesis in the early pre-rRNA processing steps at sites A0, A1 and A2 that are required for proper maturation of the 18S RNA. RCL1 has two isoforms with the molecular weight of 41 kDa and 20 kDa.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68221-1-Ig (RCL1 antibody) at dilution of 1:20000 and incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.