

À des fins de recherche uniquement

# Anticorps Monoclonal anti-COASY

Numéro de catalogue: 68358-1-Ig



## Informations de base

Numéro de catalogue:	68358-1-Ig	Numéro d'acquisition GenBank:	BC020985	Méthode de purification:	Purification par protéine G
Taille:	150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	80347	CloneNo.:	1D1C11
Hôte:	Mouse	Nom complet:	Coenzyme A synthase	Dilutions recommandées:	WB 1:5000-1:50000
Isotype:	IgG1	MW calculé:	564 aa, 62 kDa		
Immunogen Catalog Number:	AG33142	MW observés:	62 kDa		

## Applications

Applications testées:  
WB, ELISA

Spécificité de l'espèce:  
Humain, rat, souris

Contrôles positifs:

WB : cellules LNCaP, cellules HEK-293, cellules HeLa, cellules HepG2, cellules HSC-T6, cellules Jurkat, cellules K-562, cellules NIH/3T3

## Informations générales

Bifunctional coenzyme A synthase (COASY) is a bifunctional enzyme that catalyzes the fourth and fifth sequential steps of the CoA biosynthetic pathway. The fourth reaction is catalyzed by the phosphopantetheine adenylyltransferase, coded by the coaD domain; the fifth reaction is catalyzed by the dephospho-CoA kinase, coded by the coaE domain. May act as a point of CoA biosynthesis regulation. It was reported that having a key role in erythroid maturation and identified upstream and downstream metabolites of COASY as a potential treatment for anemia in patients with MDS-RS (PMID: 36857430). COASY has 2 isoforms with molecular weights of 62 and 65 kDa.

## Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.  
Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68358-1-Ig (COASY antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.