

À des fins de recherche uniquement

Anticorps Monoclonal anti-TBCC

Numéro de catalogue: 68371-1-Ig



Informations de base

Numéro de catalogue:	68371-1-Ig	Numéro d'acquisition GenBank:	BC020170	Méthode de purification:	Purification par protéine G
Taille:	150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	6903	CloneNo.:	2B7D9
Hôte:	Mouse	Nom complet:	tubulin folding cofactor C	Dilutions recommandées:	WB 1:5000-1:50000
Isotype:	IgG1	MW calculé	346 aa, 39 kDa		
Immunogen Catalog Number:	AG8962	MW observés:	39 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : HT-1376 cells, cellules 4T1, cellules HeLa, cellules HepG2, cellules HSC-T6, cellules MDA-MB-231, cellules MDA-MB-468, cellules NIH/3T3
Spécificité de l'espèce:	Humain, rat, souris		

Informations générales

TBCC, also named CFC. Expressed in the retina. Expressed in the rod and cone photoreceptors, extending from the inner segments (IS), through the outer nuclear layer (ONL) and into the synapses in the outer plexiform layer (OPL). Strongly expressed to the photoreceptor connecting cilium at the tips of the IS (at protein level). Cofactor C was cytoplasmic in distribution (PMID: 12417528).

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68371-1-Ig (TBCC antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Lamin B1 Monoclonal antibody (HRP-66095) as loading control.