

À des fins de recherche uniquement

Anticorps Monoclonal anti-ARL3

Numéro de catalogue: 68421-1-Ig



Informations de base

Numéro de catalogue:	68421-1-Ig	Numéro d'acquisition GenBank:	BC009841	Méthode de purification:	Purification par protéine A
Taille:	150ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	403	CloneNo.:	2E5B12
Hôte:	Mouse	Nom complet:	ADP-ribosylation factor-like 3	Dilutions recommandées:	WB 1:2000-1:10000
Isotype:	IgG2b	MW calculé	20 kDa		
Immunogen Catalog Number:	AG21387	MW observés:	20-24 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : cellules HeLa, tissu rétinien de rat, tissu testiculaire de lapin
Spécificité de l'espèce:	Humain, rat		

Informations générales

ARL3 (ADP-ribosylation factor-like 3) is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity. ARL3 is believed to be involved in ciliary and microtubule-dependent processes.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées



HeLa cells were subjected to SDS PAGE followed by western blot with 68421-1-Ig (ARL3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



rabbit testis tissue were subjected to SDS PAGE followed by western blot with 68421-1-Ig (ARL3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



rat retina tissue were subjected to SDS PAGE followed by western blot with 68421-1-Ig (ARL3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.