

À des fins de recherche uniquement

Anticorps Monoclonal anti-ATXN10

Numéro de catalogue: 68455-1-Ig



Informations de base

Numéro de catalogue: 68455-1-Ig	Numéro d'acquisition GenBank: BC007508	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 25814	CloneNo.: 1B4E4
Hôte: Mouse	Nom complet: ataxin 10	Dilutions recommandées: WB 1:5000-1:50000
Isotype: IgG2a	MW calculé: 53 kDa	
Immunogen Catalog Number: AG8254	MW observés: 45 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules HeLa, cellules HepG2, cellules Jurkat
Spécificité de l'espèce: Humain	

Informations générales

ATXN10, also named as SCA10, is a 475 amino acid protein, which belongs to the ataxin-10 family. ATXN10 is expressed in the central nervous system. ATXN10 is necessary for the survival of cerebellar neurons. ATXN10 may play a role in the maintenance of a critical intracellular glycosylation level and homeostasis.

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 68455-1-Ig (ATXN10 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.