

À des fins de recherche uniquement

# Anticorps Recombinant de lapin anti-Cystatin C



Numéro de catalogue: 82441-1-RR

## Informations de base

Numéro de catalogue: 82441-1-RR	Numéro d'acquisition GenBank: BC013083	Méthode de purification: Purification par protéine A
Taille: 100ul, Concentration: 750 µg/ml by Nanodrop;	Identification du gène (NCBI): 1471	CloneNo.: 1I20
Hôte: Lapin	Nom complet: cystatin C	Dilutions recommandées: WB 1:2000-1:15000 IHC 1:150-1:600
Isotype: IgG	MW calculé 146 aa, 16 kDa	
Immunogen Catalog Number: AG2890	MW observés: 13 kDa	

## Applications

Applications testées:  
IHC, WB, ELISA

Spécificité de l'espèce:  
Humain

Contrôles positifs:

WB : cellules HepG2, lait humain, urine humaine  
IHC : tissu rénal humain,

**Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (\*) À défaut, 'le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.**

## Informations générales

### Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 82441-1-RR (Cystatin C antibody) at dilution of 1:3300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 82441-1-RR (Cystatin C antibody) at dilution of 1:300 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 82441-1-RR (Cystatin C antibody) at dilution of 1:300 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).