

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG2a Isotype Control



Numéro de catalogue: **AB-65244**

Informations de base

Numéro de catalogue:

AB-65244

Taille:

100ug, 0.2 mg/ml

Hôte:

Mouse

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

MOPC-173

Excitation/Emission maxima
wavelengths:

404 nm / 458 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1×10^6 human PBMCs were surface stained with 0.2 μ g Atlantic Blue™ Anti-Mouse IgG2a Isotype Control (AB-65244, Clone:MOPC-173) (red), Atlantic Blue™ Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.