

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD33

Numéro de catalogue: **AB-65272**



Informations de base

Numéro de catalogue: AB-65272	Numéro d'acquisition GenBank: BC028152	Méthode de purification: Purification par affinité
Taille: 100tests , 5 ul/test	Identification du gène (NCBI): 945	CloneNo.: WM53
Hôte: Mouse	Nom complet: CD33 molecule	Excitation/Emission maxima wavelengths: 404 nm / 458 nm
Isotype: IgG1, kappa	MW calculé 364 aa, 40 kDa	

Applications

Applications testées:
FC
Spécificité de l'espèce:
Humain

Informations générales

CD33, also referred to as Siglec-3, is a 67-kDa glycosylated transmembrane protein of sialic acid-binding immunoglobulin-like lectin (siglec) family. CD33 is expressed on monocytes, myeloid progenitors, granulocytes, mast cells, some T cells. It may mediate cell-to-cell adhesion and act as a receptor that inhibits the proliferation of normal and leukemic myeloid cells.

Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD33 (AB-65272, Clone: WM53) (red) or unstained. Cells were not fixed. Monocytes were gated.