

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD7

Numéro de catalogue: **CL405-65203**



Informations de base

Numéro de catalogue:
CL405-65203

Taille:
100tests, 5 µl/test

Hôte:
Mouse

Isotype:
Mouse IgG2a, kappa

Numéro d'acquisition GenBank:
BC009293

Identification du gène (NCBI):
924

Nom complet:
CD7 molecule

MW calculé
240 aa, 25 kDa

Méthode de purification:
Purification par affinité

CloneNo.:
4H9

Excitation/Emission maxima
wavelengths:
399 nm / 422 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
Humain

Informations générales


Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD7 (CL405-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD7 (CL405-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.