

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TERF2IP



Numéro de catalogue: **CL488-14595**

Informations de base

Numéro de catalogue: CL488-14595	Numéro d'acquisition GenBank: BC005841	Méthode de purification: Purifié par affinité contre l'antigène
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 54386	Dilutions recommandées: IF 1:50-1:500
Hôte: Lapin	Nom complet: telomeric repeat binding factor 2, interacting protein	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG	MW calculé 44 kDa	
Immunogen Catalog Number: AG6130	MW observés: 60 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules HeLa,
Spécificité de l'espèce: Humain	

Informations générales

TERF2IP, also named as DRIP5, RAP1, is a 399 amino acid protein, which belongs to the RAP1 family. TERF2IP associates with chromosomes, both at telomeres and in extratelomeric sites. TERF2IP acts both as a regulator of telomere function and as a transcription regulator. TERF2IP is involved in the regulation of telomere length and protection as a component of the shelterin complex (telosome). In contrast to other components of the shelterin complex, it is dispensable for telomere capping and does not participate in the protection of telomeres against non-homologous end-joining (NHEJ)-mediated repair. The calculated molecular weight of TERF2IP is 44 kDa, but modified TERF2IP, is about 60 kDa. (PMID: 15838871)

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 TERF2IP antibody (CL488-14595) at dilution of 1:200, CL594-phalloidin (red).