

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HNRNPU



Numéro de catalogue: **CL488-14599**

Informations de base

Numéro de catalogue: CL488-14599	Numéro d'acquisition GenBank: BC003621	Méthode de purification: Purification par affinité contre l'antigène
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3192	Dilutions recommandées: IF 1:50-1:500
Hôte: Lapin	Nom complet: heterogeneous nuclear ribonucleoprotein U (scaffold attachment factor A)	Excitation/Emission maximum wavelengths: 493 nm / 522 nm
Isotype: IgG	MW calculé 825 aa, 91 kDa	
Immunogen Catalog Number: AG6141	MW observés: 120 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules HeLa,
Spécificité de l'espèce: Humain, rat, souris	

Informations générales

HNRNPU, also named as HNRPU, SAFA, and U21.1, is a 825 amino acid protein, which is extensively phosphorylated. HNRNPU as a component of ribonucleosomes localized in cytoplasmic mRNP granules containing untranslated mRNAs. HNRNPU is a component of the CRD-mediated complex that promotes MYC mRNA stabilization. HNRNPU Plays a role in the circadian regulation of the core clock component ARNTL/BMAL1 transcription. The calculated molecular weight of HNRNPU is 90 kDa, but the 120 kDa HNRNPU protein is an abundant nucleoplasmic phosphoprotein that can be crosslinked to pre-mRNA in intact cells by UV light (PMID: 1628625).

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 HNRNPU antibody (CL488-14599) at dilution of 1:200, CL594-phalloidin (red).