

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD11c/Integrin Alpha X



Numéro de catalogue: CL488-60258

1 Publications

Informations de base

Numéro de catalogue: CL488-60258	Numéro d'acquisition GenBank: BC038237	Méthode de purification: Purification par protéine A
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3687	CloneNo.: 2F1C10
Hôte: Mouse	Nom complet: integrin, alpha X (complement component 3 receptor 4 subunit)	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG2a	MW calculé 1169 aa, 129 kDa	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG11350		

Applications

Applications testées: IF	Contrôles positifs: IF : tissu d'amygdalite humain,
Demandes citées: IF	
Spécificité de l'espèce: Humain	
Espèces citées: souris	

Informations générales

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits. CD11c/Integrin alpha X is a 145-150 kDa type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, NK cells and dendritic cells. Integrin alpha X/beta 2 acts a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis. And this antibody is CL488(Ex/Em 488 nm/515 nm) conjugated.

Publications notables

Autrice	Pubmed ID	Journal	Application
Kirsten M Reeves	34615724	Clin Cancer Res	IF

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliqotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 CD11c/Integrin Alpha X antibody (CL488-60258, Clone: 2F1C10) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 CD11c/Integrin Alpha X antibody (CL488-60258, Clone: 2F1C10) at dilution of 1:200.