

À des fins de recherche uniquement

Anticorps Monoclonal anti-IL-28A

Numéro de catalogue: CL488-60270



Informations de base

Numéro de catalogue: CL488-60270	Numéro d'acquisition GenBank: BC113583	Méthode de purification: Purification par protéine G
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 282616	CloneNo.: 3B4A1
Hôte: Mouse	Nom complet: interleukin 28A (interferon, lambda 2)IF 1:300-1:1200	Dilutions recommandées:
Isotype: IgG1	MW calculé 200 aa, 22 kDa	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG18525		

Applications

Applications testées: IF	Contrôles positifs: IF : tissu d'amygdales humain,
Spécificité de l'espèce: Humain	

Informations générales

IL28A, also named as IFNL2 and ZCYT020, is a cytokine with immunomodulatory activity. It up-regulates MHC class I antigen expression. IL28A is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IFNLR1. The ligand/receptor complex seems to signal through the Jak-STAT pathway. It plays a significant role in the antiviral immune defense in the intestinal epithelium.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-28A antibody (CL488-60270, Clone: 3B4A1) at dilution of 1:600.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-28A antibody (CL488-60270, Clone: 3B4A1) at dilution of 1:600.