

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD7

Numéro de catalogue: **CL488-65203**



## Informations de base

Numéro de catalogue:

CL488-65203

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC009293

Identification du gène (NCBI):

924

Nom complet:

CD7 molécule

MW calculé

240 aa, 25 kDa

Méthode de purification:

Purification par protéine A

CloneNo.:

4H9

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



1X10<sup>6</sup> human PBMCs were surface co-stained with APC Anti-Human CD4 and 5 ul CoraLite® Plus 488 Anti-Human CD7 (CL488-65203, Clone:4H9) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD7 (CL488-65203, Clone:4H9) (red) or Mouse IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.