

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG1 Isotype Control



Numéro de catalogue: **CL488-65212**

Informations de base

Numéro de catalogue:

CL488-65212

Taille:

100ug, 0.5 mg/mL

Hôte:

Rat

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

HRPN

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

The HRPN immunoglobulin is useful as an isotype-matched control. The HRPN immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG1 antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Mouse CD4 (CL647-65104, Clone: GK1.5) and 5 ul CoraLite® Plus 488-conjugated Anti-Mouse CD25 (CL488-65137, Clone: PC61.5). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Mouse CD4 (CL647-65104, Clone: GK1.5) and CoraLite® Plus 488-conjugated Rat IgG1 Isotype Control (CL488-65212, Clone: HRPN). Cells were not fixed.