

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD43

Numéro de catalogue: **CL488-65237**



## Informations de base

Numéro de catalogue:

CL488-65237

Taille:

100ug, 0.5 mg/ml

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

24796

Nom complet:

sialophorin

MW calculé

38 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

W3/13

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

rat

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 488 Anti-Rat CD43 (CL488-65237, Clone:W3/13), or 0.5 ug Control Antibody. Cells were not fixed.