

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD96

Numéro de catalogue: **CL488-65250**



Informations de base

Numéro de catalogue:

CL488-65250

Taille:

100tests , 5 ul/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC020749

Identification du gène (NCBI):

10225

Nom complet:

CD96 molecule

MW calculé

402 aa, 46 kDa

Méthode de purification:

N/A

CloneNo.:

NK92.39

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface co-stained with CoraLite® Plus 647 Anti-Human CD19 and 5 ul CoraLite® Plus 488 Anti-Human CD96 (CL488-65250, Clone: NK92.39) or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD96 (CL488-65250, Clone: NK92.39) (red) or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.