

À des fins de recherche uniquement

Anticorps Monoclonal anti-GNAS

Numéro de catalogue: CL488-66253



Informations de base

| | | |
|---|---|---|
| Numéro de catalogue: CL488-66253 | Numéro d'acquisition GenBank: BC002722 | Méthode de purification: Purification par protéine A |
| Taille: 100ul, Concentration: 1000 µg/ml by Nanodrop; | Identification du gène (NCBI): 2778 | CloneNo.: 1A9C10 |
| Hôte: Mouse | Nom complet: GNAS complex locus | Dilutions recommandées: IF 1:50-1:500 |
| Isotype: IgG2b | MW calculé 45 kDa | Excitation/Emission maxima wavelengths: 493 nm / 522 nm |
| Immunogen Catalog Number: AG20530 | MW observés: 46 kDa | |

Applications

| | |
|---|---|
| Applications testées: FC (Intra), IF | Contrôles positifs: IF : cellules MCF-7, |
| Spécificité de l'espèce: Humain | |

Informations générales

Guanine nucleotide binding protein (G protein), alpha stimulating activity polypeptide 1 (GNAS1) is the ubiquitously expressed heterotrimeric G protein that couples receptors to the effector enzyme adenylyl cyclase and is required for receptor-stimulated intracellular cAMP generation. Mutations of Gs(alpha) residues involved in the GTPase reaction that lead to constitutive activation are present in endocrine tumors, fibrous dysplasia of bone, and McCune-Albright syndrome. The molecular weight of Gs(alpha) protein is about 46 kDa.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using CoraLite® Plus 488 GNAS antibody (CL488-66253, Clone: 1A9C10) at dilution of 1:200.



1X10⁶ MCF-7 cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human GNAS (CL488-66253, Clone:1A9C10) (red), or 0.4 ug Mouse IgG2b Isotype Control (CL488-66360-3, Clone: K11B8C4B5) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).