

À des fins de recherche uniquement

# Anticorps Monoclonal anti-SMARCA4

Numéro de catalogue: CL488-66561



## Informations de base

<b>Numéro de catalogue:</b> CL488-66561	<b>Numéro d'acquisition GenBank:</b> BC150298	<b>Méthode de purification:</b> Purification par protéine A
<b>Taille:</b> 100ul, Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	<b>Identification du gène (NCBI):</b> 6597	<b>CloneNo.:</b> 2E6B6
<b>Hôte:</b> Mouse	<b>Nom complet:</b> SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4	<b>Dilutions recommandées:</b> IF 1:50-1:500
<b>Isotype:</b> IgG2b	<b>MW calculé</b> 1647 aa, 185 kDa	<b>Excitation/Emission maxima wavelenghts:</b> 493 nm / 522 nm
<b>Immunogen Catalog Number:</b> AG16256		

## Applications

<b>Applications testées:</b> IF	<b>Contrôles positifs:</b> IF : cellules HepG2,
<b>Spécificité de l'espèce:</b> Humain	

## Informations générales

SMARCA4, also named as BAF190A, BRG1, SNF2B and SNF2L4, belongs to the SNF2/RAD54 helicase family. SMARCA4 is a transcriptional coactivator cooperating with nuclear hormone receptors to potentiate transcriptional activation. It is a component of the CREST-BRG1 complex, a multiprotein complex that regulates promoter activation by orchestrating a calcium-dependent release of a repressor complex and a recruitment of an activator complex. It is also involved in vitamin D-coupled transcription regulation via its association with the WINAC complex, a chromatin-remodeling complex recruited by vitamin D receptor (VDR), which is required for the ligand-bound VDR-mediated transrepression of the CYP27B1 gene. And this antibody is CL488 Ex/Em 488 nm/515 nm conjugated.

## Stockage

**Stockage:**  
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.  
**Tampon de stockage:**  
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CL488-66561 (SMARCA4 antibody) at dilution of 1:50.