

À des fins de recherche uniquement

# Anticorps Monoclonal anti-HDAC9

Numéro de catalogue: **CL488-67364**



## Informations de base

Numéro de catalogue:	CL488-67364	Numéro d'acquisition GenBank:	BC152405	Méthode de purification:	Purification par protéine A
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	9734	CloneNo.:	1G6C5
Hôte:	Mouse	Nom complet:	histone deacetylase 9	Excitation/Emission maxima wavelengths:	493 nm / 522 nm
Isotype:	IgG2a	MW calculé	1011 aa, 111 kDa		
Immunogen Catalog Number:	AG28514	MW observés:	130 kDa		

## Applications

Applications testées:  
FC (Intra)  
Spécificité de l'espèce:  
Humain

## Informations générales

HDAC9, also named as Histone deacetylase 7B, is a 1011 amino acid protein, which is responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. HDAC9 represses MEF2-dependent transcription. HDAC9 is broadly expressed, with highest levels in brain, heart, muscle and testis.

## Stockage

Stockage:  
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.  
Tampon de stockage:  
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



1X10<sup>6</sup> HeLa cells were intracellularly stained with 0.4 ug Coralite® Plus 488 Anti-Human HDAC9 (CL488-67364, Clone:1G6C5) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).