

À des fins de recherche uniquement

Anticorps Monoclonal anti-VASH2

Numéro de catalogue: **CL488-67753**



Informations de base

Numéro de catalogue:	CL488-67753	Numéro d'acquisition GenBank:	BC028194	Méthode de purification:	Purification par protéine G
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	79805	CloneNo.:	3A11D10
Hôte:	Mouse	Nom complet:	vasohibin 2	Dilutions recommandées:	IF 1:50-1:500
Isotype:	IgG1	MW calculé	355 aa, 40 kDa	Excitation/Emission maxima wavelengths:	493 nm / 522 nm
Immunogen Catalog Number:	AG29985	MW observés:	34 kDa		

Applications

Applications testées:	IF	Contrôles positifs:	IF : cellules SKOV-3,
Spécificité de l'espèce:	Humain		

Informations générales

VASH2, also named as vasohibin 2, is an important pro-angiogenesis factor in solid tumor. It has been reported that VASH2 is expressed in mononuclear cells mobilized from bone marrow to promote angiogenesis. VASH2 also plays a key role in axon formation. VASH2 is localized as a nuclear type (with 311 amino acid residues) and cytoplasmic type (with 355 amino acid residues and low abundance). Cytoplasmic VASH2 is associated with carcinoma angiogenesis, while nuclear VASH2 may be associated with cell proliferation. 67753-1-Ig antibody recognizes both the nuclear and cytoplasmic isoforms. (PMID:23615928; 19204325; 31235911; 26177649)

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed SKOV-3 cells using CoraLite® Plus 488 VASH2 antibody (CL488-67753, Clone: 3A11D10) at dilution of 1:200.