

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-C1QBP



Numéro de catalogue: **CL555-24474** **Phare**

Informations de base

Numéro de catalogue: CL555-24474	Numéro d'acquisition GenBank: BC013731	Méthode de purification: Purification par affinité contre l'antigène
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 708	Dilutions recommandées: IF 1:50-1:500
Hôte: Lapin	Nom complet: complement component 1, q subcomponent binding protein	Excitation/Emission maxima wavelengths: 557 nm / 570nm
Isotype: IgG	MW calculé: 282 aa, 31 kDa	
Immunogen Catalog Number: AG19773	MW observés: 32 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules NIH/3T3,
Spécificité de l'espèce: Humain, rat, souris	

Informations générales

C1QBP, also named as gC1q receptor (gC1qR), p32, p33, and hyaluronan-binding protein 1 (HABP1), is a protein initially copurified with splicing factor SF2 (PMID: 1830244). The protein is synthesized as a pro-protein of 282 amino acids (aa) that is post-translationally processed by removal of the initial 73 aa to a mature protein of 209 aa (PMID: 8262387). C1QBP is an evolutionary conserved and ubiquitously expressed multifunctional protein and has been reported to be a predominantly mitochondrial matrix protein involved in inflammation and infection processes, mitochondrial ribosome biogenesis, regulation of apoptosis and nuclear transcription, and pre-mRNA splicing (PMID: 28942965).

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed NIH/3T3 cells using CoraLite®555 C1QBP antibody (CL555-24474) at dilution of 1:100.