

À des fins de recherche uniquement

Anticorps Monoclonal anti-pan Ras

Numéro de catalogue: **CL555-60309** **Phare**



Informations de base

Numéro de catalogue: CL555-60309	Numéro d'acquisition GenBank: BC013572	Méthode de purification: Purification par protéine G
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3845	CloneNo.: 4H4G7
Hôte: Mouse	Nom complet: v-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog	Excitation/Emission maxima wavelengths: 557 nm / 570nm
Isotype: IgG1	MW calculé: 188 aa, 21 kDa	
Immunogen Catalog Number: AG2700	MW observés: 21 kDa	

Applications

Applications testées:
FC (Intra)

Spécificité de l'espèce:
Humain, porc, rat, souris

Informations générales

The 21 kDa guanine-nucleotide binding proteins (K-Ras, H-Ras, and N-Ras) belong to the Ras oncogene family, whose members are related to the transforming genes of mammalian sarcoma retroviruses. K-Ras, H-Ras, and N-Ras have similar structure and sequences. These proteins can bind GTP and GDP, and they have intrinsic GTPase activity. The ras genes are ubiquitously expressed although mRNA analysis suggests different level expression in tissue. Mutations in each ras gene frequently were found in different tumors, suggesting their involvement in the development of specific neoplasia. This antibody can recognize K-Ras, H-Ras, and N-Ras.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.


Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ HeLa cells were intracellularly stained with 0.4 ug Coralite®555 Anti-Human pan Ras (CL555-60309, Clone:4H4G7) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).