

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG2a Isotype Control



Numéro de catalogue: **CL555-65209**

Informations de base

Numéro de catalogue:
CL555-65209

Taille:
100ug, 0.5 mg/ml

Hôte:
Rat

Isotype:
IgG2a, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

2A3

Excitation/Emission maxima
wavelengths:

554 nm / 570 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1×10^6 mouse splenocytes were surface stained with 0.5 μ g CoraLite® Plus 555 Rat IgG2a Isotype Control (CL555-65209, Clone: 2A3) (red), or 0.5 μ g Control Antibody (green), or unstained (blue). Cells were not fixed.