

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD48

Numéro de catalogue: **CL555-65227**



## Informations de base

Numéro de catalogue:  
CL555-65227

Taille:  
100ug, 0.5 mg/ml

Hôte:  
Armenian Hamster

Isotype:  
IgG

Numéro d'acquisition GenBank:  
BC060977

Identification du gène (NCBI):  
12506

Nom complet:  
CD48 antigen

Méthode de purification:  
Purification par affinité

CloneNo.:  
HM48-1

Excitation/Emission maxima  
wavelengths:  
554 nm / 570 nm

## Applications

Applications testées:  
FC

Spécificité de l'espèce:  
souris

## Informations générales

### Stockage

Stockage:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Tampon de stockage:  
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug CoraLite® Plus 555 Anti-Mouse CD48 (CL555-65227, Clone:HM48-1) (red), or 0.5 ug Isotype Control. Cells were not fixed.