

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-HuR



Numéro de catalogue: **CL594-11910**

Phare

## Informations de base

Numéro de catalogue:

CL594-11910

Taille:

100ul , Concentration: 1000 µg/ml by 1994 Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG2249

Numéro d'acquisition GenBank:

BC003376

Identification du gène (NCBI):

1994

Nom complet:

ELAV (embryonic lethal, abnormal vision, Drosophila)-like 1 (Hu antigen R)

MW calculé

326 aa, 36 kDa

MW observés:

30-37 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

## Applications

Applications testées:

FC (Intra)

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

ELAVL1, also named as HUR, belongs to the RRM elav family. It is involved in 3'-UTR ARE-mediated MYC stabilization. ELAVL1 binds avidly to the AU-rich element in FOS and IL3/interleukin-3 mRNAs. In the case of the FOS AU-rich element, ELAVL1 binds to a core element of 27 nucleotides that contain AUUUA, AUUUUA and AUUUUUA motifs.

## Stockage

Stockage:

Stocker à -20 °C. Éviter toute exposition à la lumière.

Tampon de stockage:

PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> MCF-7 cells were intracellularly stained with 0.4 ug CoraLite® 594 Anti-Human HuR (CL594-11910) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).