

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-TBR1



Numéro de catalogue: **CL594-20932**

Informations de base

Numéro de catalogue:	CL594-20932	Numéro d'acquisition GenBank:	BC104844	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	10716	Dilutions recommandées:	IF 1:50-1:500
Hôte:	Lapin	Nom complet:	T-box, brain, 1	Excitation/Emission maxima wavelengths:	588 nm / 604 nm
Isotype:	IgG	MW calculé	682 aa, 74 kDa		
Immunogen Catalog Number:	AG14935	MW observés:	74 kDa		

Applications

Applications testées:	IF	Contrôles positifs:	IF : tissu cérébral de souris,
Spécificité de l'espèce:	Humain, rat, souris		

Informations générales

TBR1, also named as T-box brain protein 1, is a 682 amino acid protein, which contains one T-box DNA-binding domain and localizes in the nucleus. TBR1 is expressed in the brain and as a transcriptional regulator is involved in developmental processes. TBR1 is required for normal brain development.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CoraLite® 594 TBR1 antibody (CL594-20932) at dilution of 1:200.