

À des fins de recherche uniquement

Anticorps Monoclonal anti-Gastrin

Numéro de catalogue: **CL594-60346**



Informations de base

Numéro de catalogue:
CL594-60346

Taille:
100ul , Concentration: 1000 µg/ml by
Nanodrop;

Hôte:
Mouse

Isotype:
IgG2a

Immunogen Catalog Number:
AG12795

Numéro d'acquisition GenBank:
BC069724

Identification du gène (NCBI):
2520

Nom complet:
gastrin

MW calculé
101 aa, 11 kDa

Méthode de purification:
Purification par protéine A

CloneNo.:
7G9C3

Excitation/Emission maxima
wavelengths:
588 nm / 604 nm

Applications

Applications testées:
Spécificité de l'espèce:
Humain

Informations générales

Gastrin is a hormone produced primarily by G-cells in the stomach, where it functions to stimulate acid secretion by gastric parietal cells. Gastrin is a promiscuous hormone. Thus, the gastrin gene is expressed in common cancers such as bronchogenic, colorectal, ovarian and pancreatic carcinomas.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

