

À des fins de recherche uniquement

Anticorps Monoclonal anti-KLRG1

Numéro de catalogue: **CL594-65084**



Informations de base

Numéro de catalogue:
CL594-65084

Taille:
100ug, 0.5 mg/ml

Hôte:
Golden Syrian Hamster

Isotype:
IgG

Numéro d'acquisition GenBank:
BC101953

Identification du gène (NCBI):
50928

Nom complet:
killer cell lectin-like receptor
subfamily G, member 1

Méthode de purification:
Purification par affinité

CloneNo.:
2F1

Excitation/Emission maxima
wavelengths:
588 nm / 604 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
souris

Informations générales

KLRG1 (killer cell lectin-like receptor subfamily G, member 1) is a C-type lectin inhibitory receptor that contains an immunoreceptor tyrosine-based inhibitory motif (ITIM) in its cytoplasmic domain. KLRG1 is expressed by subsets of NK and T cells, existing both as a monomer and as a disulfide-linked homodimer (PMID: 10741410; 19604491). It is considered to be a cell differentiation marker for NK and T cells and is strongly induced by viral and other infections (PMID: 20373518). Through interactions with members of the cadherin family, KLRG1 plays an inhibitory role on NK- and T-cell function (PMID: 19009530; 20373518).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.5 ug CoraLite®594 Anti-Mouse KLRG1 (CL594-65084, Clone: 2F1), or 0.5 ug Isotype Control. Cells were not fixed.