

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD48

Numéro de catalogue: **CL594-65227**



Informations de base

Numéro de catalogue:
CL594-65227

Taille:
100ug, 0.5 mg/ml

Hôte:
Armenian Hamster

Isotype:
IgG

Numéro d'acquisition GenBank:
BC060977

Identification du gène (NCBI):
12506

Nom complet:
CD48 antigen

Méthode de purification:
Purification par protéine A

CloneNo.:
HM48-1

Excitation/Emission maxima
wavelengths:
588 nm / 604 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
souris

Informations générales


Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.5 ug Coralite® 594 Anti-Mouse CD48 (CL594-65227, Clone:HM48-1) (red), or 0.5 ug Isotype Control. Cells were not fixed.