

À des fins de recherche uniquement

Anticorps Monoclonal anti-HSP60

Numéro de catalogue: CL594-66041

1 Publications



Informations de base

Numéro de catalogue: CL594-66041	Numéro d'acquisition GenBank: BC003030	Méthode de purification: Purification par protéine G
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3329	CloneNo.: 2F10E7
Hôte: Mouse	Nom complet: heat shock 60kDa protein 1 (chaperonin)	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG1	MW calculé: 61 kDa	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG7848	MW observés: 60 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules HepG2,
Demandes citées: FC	
Spécificité de l'espèce: Humain, rat, souris	
Espèces citées: souris	

Informations générales

HSPD1, also known as HSP60, belongs to the chaperonin family and acts as a chaperone to enhance cell survival under physiological stresses. Hsp60 has been shown to be connected with many aspects of cell functions such as protein folding and assembling of polypeptide chains in mitochondria. Recently it has been reported that HSP60 is associated with apoptosis or inhibition of cancer cell growth. (21822415)

Publications notables

Autrice	Pubmed ID	Journal	Application
Yavuz F Yazicioglu	37095377	Nat Immunol	FC

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using CL594-66041 (HSP60 antibody) at dilution of 1:50.