

À des fins de recherche uniquement

Anticorps Monoclonal anti-HE4

Numéro de catalogue: CL594-66557



Informations de base

Numéro de catalogue: CL594-66557	Numéro d'acquisition GenBank: BC046106	Méthode de purification: Purification par protéine G
Taille: 100ul, Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 10406	CloneNo.: 5B2D9
Hôte: Mouse	Nom complet: WAP four-disulfide core domain 2	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG1	MW calculé 13 kDa	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG5936		

Applications

Applications testées: IF	Contrôles positifs: IF : cellules SKOV-3,
Spécificité de l'espèce: Humain	

Informations générales

WFDC2, also named as HE4, is a member of the WFDC domain family. The WFDC domain, or WAP Signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many family members. HE4 is expressed in a number of normal tissues, and it is also highly expressed in a number of tumors cells lines, such ovarian, colon, breast, lung and renal cells lines. In some study, HE4 was referred to as a biomarker for ovarian carcinoma.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed SKOV-3 cells using CoraLite® 594 HE4 antibody (CL594-66557, Clone: 5B2D9) at dilution of 1:200.