

À des fins de recherche uniquement

Anticorps Monoclonal anti-PARK2/Parkin



Numéro de catalogue: **CL594-66674**

Informations de base

Numéro de catalogue:

CL594-66674

Taille:

100ul

Hôte:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG5179

Numéro d'acquisition GenBank:

BC022014

Identification du gène (NCBI):

5071

Nom complet:

Parkinson disease (autosomal recessive, juvenile) 2, parkin

MW calculé

52 kDa

MW observés:

42-52 kDa

Méthode de purification:

Purification par protéine A

CloneNo.:

2H5A7

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Applications testées:

Spécificité de l'espèce:

Humain, souris

Informations générales

Parkin, a RING-type E3 ubiquitin-protein ligase, is involved in the ubiquitination pathway and contributes to protection from neurotoxicity induced by unfolded protein stresses. Its ubiquitin-protein ligase activity promotes the degradation of a variety of proteins including itself. Mutations in Parkin are implicated in the pathogenesis of autosomal recessive familial Parkinson's disease. It has 8 isoforms produced by alternative splicing.

Stockage

Stockage:

Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

