

À des fins de recherche uniquement

# Anticorps Monoclonal anti-NRBP2

Numéro de catalogue: **CL594-66687**



## Informations de base

|                           |                                                |                                |                                    |                                         |                             |
|---------------------------|------------------------------------------------|--------------------------------|------------------------------------|-----------------------------------------|-----------------------------|
| Numéro de catalogue:      | CL594-66687                                    | Numéro d'acquisition GenBank:  | BC113873                           | Méthode de purification:                | Purification par protéine A |
| Taille:                   | 100ul , Concentration: 1000 µg/ml by Nanodrop; | Identification du gène (NCBI): | 340371                             | CloneNo.:                               | 4F10H7                      |
| Hôte:                     | Mouse                                          | Nom complet:                   | nuclear receptor binding protein 2 | Dilutions recommandées:                 | IF 1:50-1:500               |
| Isotype:                  | IgG2a                                          | MW calculé                     | 501 aa, 58 kDa                     | Excitation/Emission maxima wavelengths: | 588 nm / 604 nm             |
| Immunogen Catalog Number: | AG17011                                        | MW observés:                   | 58 kDa, 33 kDa                     |                                         |                             |

## Applications

|                           |                     |
|---------------------------|---------------------|
| Applications testées:     | Contrôles positifs: |
| IF                        | IF : cellules HeLa, |
| Spécificité de l'espèce:  |                     |
| Humain, porc, rat, souris |                     |

## Informations générales

Nuclear receptor binding protein 2 (NRBP2), also named as TRG16, belongs to the Ser/Thr protein kinase family. NRBP2 can regulate apoptosis of neural progenitor cells during their differentiation. Expression of NRBP2 was up-regulated following differentiation in rodent neural stem/progenitor cells. It has 4 isoforms with MW of 58, 33, 36 and 14 kDa. Endogenous NRBP2 in HEK293 cells has an apparent molecular mass of 55-60 kD by SDS-PAGE.

## Stockage

**Stockage:**  
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.  
**Tampon de stockage:**  
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.  
L'aliquotage n'est pas nécessaire pour le stockage à -20C


\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées

---



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CL594-66687 (NRBP2 antibody) at dilution of 1:100.