

À des fins de recherche uniquement

Anticorps Monoclonal anti-HSPH1

Numéro de catalogue: **CL594-66723**



Informations de base

Numéro de catalogue:	CL594-66723	Numéro d'acquisition GenBank:	BC037553	Méthode de purification:	Purification par protéine G
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	10808	CloneNo.:	3H5A5
Hôte:	Mouse	Nom complet:	heat shock 105kDa/110kDa protein 1	Dilutions recommandées:	IF 1:50-1:500
Isotype:	IgG1	MW calculé	858 aa, 97 kDa	Excitation/Emission maxima wavelengths:	588 nm / 604 nm
Immunogen Catalog Number:	AG27521	MW observés:	110 kDa		

Applications

Applications testées:	Contrôles positifs:
IF	IF : cellules HeLa,
Spécificité de l'espèce:	
Humain	

Informations générales

HSP105, also known as HSP110 or HSPH1, belongs to the heat shock protein (HSP) family. Human HSP105 is a high-molecular-weight chaperone protein expressed at constitutively low levels as a cytoplasmic α -isoform and as an inducible nuclear β -isoform on exposure to various forms of stress. HSP105 is constitutively overexpressed in several solid tumors, including melanoma, breast, thyroid, and gastroenteric cancers, and exerts antiapoptotic functions. Recently HSP105 has been identified as a novel candidate biomarker of lymphoma aggressiveness. Western blot analysis using this antibody detected the phosphorylated protein around 110 kDa in various lysates.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CL594-66723 (HSPH1 antibody) at dilution of 1:100.