

À des fins de recherche uniquement

Anticorps Monoclonal anti-PRMT2

Numéro de catalogue: CL594-66885



Informations de base

Numéro de catalogue: CL594-66885	Numéro d'acquisition GenBank: BC000727	Méthode de purification: Purification par protéine A
Taille: 100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI): 3275	CloneNo.: 1F1B9
Hôte: Mouse	Nom complet: protein arginine methyltransferase 2	Dilutions recommandées: IF 1:50-1:500
Isotype: IgG2a	MW calculé 49 kDa	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Immunogen Catalog Number: AG1806	MW observés: 45-50 kDa	

Applications

Applications testées: IF	Contrôles positifs: IF : cellules HeLa,
Spécificité de l'espèce: Humain, rat, souris	

Informations générales

PRMT2 (also known as HRMT1L1) is a member of the arginine methyltransferase family which catalyzes arginine methylation and regulates diverse cellular processes, including transcription, translation, DNA repair, and RNA processing. Recently PRMT2 has been identified as a pro-tumorigenic factor in malignant gliomas.

Stockage

Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (-20°C Methanol) fixed HeLa cells using CoraLite®594 PRMT2 antibody (CL594-66885, Clone: 1F1B9) at dilution of 1:200.