

À des fins de recherche uniquement

Anticorps Monoclonal anti-BSEP

Numéro de catalogue: **CL594-67512**



Informations de base

Numéro de catalogue:	CL594-67512	Numéro d'acquisition GenBank:	NM_003742	Méthode de purification:	Purification par protéine G
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	8647	CloneNo.:	3C11D5
Hôte:	Mouse	Nom complet:	ATP-binding cassette, sub-family B (MDR/TAP), member 11	Dilutions recommandées:	IF 1:500-1:2000
Isotype:	IgG1	MW calculé	146 kDa	Excitation/Emission maxima wavelengths:	588 nm / 604 nm
Immunogen Catalog Number:	AG29135	MW observés:	150-160 kDa		

Applications

Applications testées:	Contrôles positifs:
IF	IF : tissu hépatique de souris,
Spécificité de l'espèce:	
Humain	

Informations générales

ABCB11, also named as BSEP, belongs to the ABC transporter superfamily. ABCB11 involved in the ATP-dependent secretion of bile salts into the canaliculus of hepatocytes. Mutation of ABCB11 will cause the progressive familial intrahepatic cholestasis type 2 (PFIC2) and the benign recurrent intrahepatic cholestasis type 2 (BRIC2). Genetic variations in ABCB11 may play a role in drug-induced cholestasis causing liver damage.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed mouse liver tissue using CoraLite®594 BSEP antibody (CL594-67512, Clone: 3C11D5) at dilution of 1:1000.