

À des fins de recherche uniquement

Anticorps Monoclonal anti-HNRNPA1

Numéro de catalogue: **CL594-67844**



Informations de base

Numéro de catalogue:	CL594-67844	Numéro d'acquisition GenBank:	BC009600	Méthode de purification:	Purification par protéine A
Taille:	100ul , Concentration: 1000 µg/ml by Nanodrop;	Identification du gène (NCBI):	3178	CloneNo.:	1A5A4
Hôte:	Mouse	Nom complet:	heterogeneous nuclear ribonucleoprotein A1	Dilutions recommandées:	IF 1:400-1:1600
Isotype:	IgG2b	MW calculé	34 kDa	Excitation/Emission maxima wavelengths:	588 nm / 604 nm
Immunogen Catalog Number:	AG17415	MW observés:	34 kDa		

Applications

Applications testées:	Contrôles positifs:
IF	IF : cellules HeLa,
Spécificité de l'espèce:	
Humain, souris	

Informations générales

In eukaryotic cells, nascent RNA polymerase II transcripts are associated in the nucleus with specific proteins to form ribonucleoprotein complexes called HNRNP or 40S. Protein moiety of the HNRNP particle has 6 major components called core proteins, and HNRPA1 is one of them. hnRNP A1 is a multifunctional RNA binding protein which is involved in pre-mRNA splicing, export of mature transcripts from the nucleus, mRNA turnover, and internal ribosome entry site (IRES)-mediated translation. HNRNPA1 has three isoforms with MW 30 kDa, 34 kDa and 39 kDa.

Stockage


Stockage:
Stocker à -20 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite®594 HNRNPA1 antibody (CL594-67844, Clone: 1A5A4) at dilution of 1:800.