

À des fins de recherche uniquement

Anticorps Monoclonal anti-IL-12/IL-23 p40



Numéro de catalogue: **CL647-65097**

Informations de base

Numéro de catalogue:
CL647-65097

Taille:
100tests , 5 µl/test

Hôte:
Rat

Isotype:
Rat IgG2a, kappa

Numéro d'acquisition GenBank:
BC103608

Identification du gène (NCBI):
16160

Nom complet:
interleukin 12b

Méthode de purification:
Purification par affinité

CloneNo.:
C17.8

Excitation/Emission maxima
wavelengths:
654 nm / 674 nm

Applications

Applications testées:
FC (Intra)

Spécificité de l'espèce:
souris

Informations générales


Stockage

Stockage:
Stocker à une température comprise entre 2 °C et 8 °C. Éviter toute exposition à la lumière.
Tampon de stockage:
PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



IFN- γ primed Balb/c mouse peritoneal macrophages stimulated with LPS in the presence of a protein transport inhibitor were treated with fixation and permeabilization reagents and intracellularly stained with CoraLite® Plus 647-Rat IgG2a Isotype Control (left) or 5 μ l CoraLite® Plus 647 Anti-Mouse IL-12/IL-23 p40 (CL647-65097, Clone:C17.8) (right).